

# Chicken Anti-Goat IgG (H+L)

Chicken anti-goat IgG (H L) secondary antibody labeled with our superior CF  $\!\!\!\! ^{^{^\circ}}$  Dyes or HRP.



## **Product Description**

This is a chicken anti-goat IgG (H L) secondary antibody labeled with our superior CF® Dyes or HRP.

- Available in 7 bright and photostable CF® Dyes, and HRP
- Suitable for western, immunofluorescence, and immunohistology in FFPE tissues

**Note:** Conjugates of blue fluorescent dyes like CF®350 are not recommended for detecting low abundance targets, because blue dyes have lower fluorescence and can give higher non-specific background than other dye colors.

Call us: 800-304-5357 Email: btinfo@biotium.com

### **Product attributes**

Antibody type	Secondary		
Clonality	Polyclonal		
Host species	Chicken		
Antibody reactivity (target)	Goat IgG		
Species reactivity	Goat		
Cross adsorption	Not cross-adsorbed		
Concentration	2 mg/mL (CF® dye conjugates), 1 mg/mL (HRP, AP conjugates)		
Antibody/conjugate formulation	Liquid: PBS/50% glycerol/2 mg/mL BSA/0.05% azide, Lyophilized: PBS/15 mg/mL BSA/20 mg/mL trehalose after reconstitution, HRP conjugates: PBS/50% glycerol/15 mg/mL BSA, HRP conjugates (lyophilized): PBS/10 mg/mL trehalose after reconstitution		
Secondary/tag antibody applications	ELISA, Flow cytometry, IHC, IF (cells or tissue sections), Western blot		
Product origin	Product may contain either bovine serum albumin (BSA) from bovine serum (Bos taurus), or recombinant BSA produced in Chinese hamster ovary cells. Inquire for the specific lot.		

## Chicken Anti-Goat IgG (H+L)

Conjugation	Ex/Em	Size	Catalog No.	Dye Features
CF®350	347/448 nm	50 uL (100 ug)	20364-50uL	CF®350 Features
		0.5 mL (1 mg)	20364-500uL	
		1 mg (lyophilized)	20364-1mg	
CF®488A	490/515 nm	50 uL (100 ug)	<u>20225-1</u>	CF®488A Features
		0.5 mL (1 mg)	<u>20225</u>	
		1 mg (lyophilized)	20225-1mg	
CF®543	541/560 nm	50 uL (100 ug)	20333-50uL	CF®543 Features
		0.5 mL (1 mg)	20333-500uL	
		1 mg (lyophilized)	20333-1mg	
CF®568	562/583 nm	50 uL (100 ug)	20337-50uL	CF®568 Features
		0.5 mL (1 mg)	20337-500uL	
		1 mg (lyophilized)	20337-1mg	
CF®583R	585/609 nm	50 uL (100 ug)	20889-50uL	CF®583R Features
		0.5 mL (1 mg)	20889-500uL	
CF®594	593/614 nm	50 uL (100 ug)	<u>20226-1</u>	CF®594 Features
		0.5 mL (1 mg)	<u>20226</u>	
		1 mg (lyophilized)	20226-1mg	
CF®633	630/650 nm	50 uL (100 ug)	<u>20227-1</u>	CF®633 Features
		0.5 mL (1 mg)	<u>20227</u>	
		1 mg (lyophilized)	20227-1mg	
HRP	N/A	100 uL (100 ug)	20839-100uL	
		1 mL (1 mg)	20839-1mL	
		1 mg (lyophilized)	20839-1mg	

View our full selection of <u>Secondary Antibodies</u>, or search our catalog using our <u>Antibody Finder</u>. Alternatively, you can view our <u>secondary antibody product listings</u> with catalog numbers.

CF® Dyes offer exceptional brightness and photostability. For more information see our CF® Dye technology page.

### Storage and Handling

Liquid format: Store at -20°C, protected from light. Product is stable for at least 6 months from date of receipt when stored as recommended. Liquid format antibodies contain 50% glycerol and will not freeze at -20°C.

Lyophilized format: Store at -20°C, protected from light. Product is stable for at least 6 months from date of receipt when stored as recommended. Reconstitute antibodies in water using the indicated volumes below:

CF® Dye and biotin conjugates: add 0.5 mL dH $_2\text{O}$ 

HRP or DNP conjugates: add 1 mL dH<sub>2</sub>O

Add the indicated volume of water directly to the vial containing the lyophilized antibody and mix gently to dissolve. Store reconstituted antibody at -20°C and protect from light. Aliquot to avoid repeated freeze/thaw cycles. Alternatively, an equal volume of glycerol can be mixed with the reconstituted antibody so that it will remain liquid at -20°C.

Optional: A preservative such as 0.05% sodium azide (final concentration) can be added to CF® Dye and biotin conjugates. Do not add sodium azide to HRP conjugates.

Note: Storage of the antibody for more than a day at final working dilution is not recommended.

CF is a registered trademark of Biotium, Inc.

## References

Download a list of CF® dye references.

This datasheet was generated on April 20, 2025 at 09:16:23 AM. Visit product page to check for updated information before use. Product link: <a href="http://54.245.69.9/product/chicken-anti-goat-igg-hl/">http://54.245.69.9/product/chicken-anti-goat-igg-hl/</a>